



LCMS Troubleshooting Tips



Changes in Sensitivity

Causes	Solutions
MS	<ul style="list-style-type: none">- If all peaks have changed in sensitivity check maintenance requirements; and clean the ion source and if required ion optics.- Ensure appropriate and accurate acquisition parameters are being used.- Check for adducts such as Na⁺, NH₄⁺, K⁺. Consider different glassware, make up new mobile phases.- Increase additive or buffer concentration if required.- Ensure probe is at correct distance from orifice.- Check for matrix effects such as ion suppression. Consider different LC methods and sample preparation procedures.- Replace capillary and desolvation line.- Perform and check MS tune.
Sample	<ul style="list-style-type: none">- Check for sample degradation. Prepare a fresh sample.- Check sample preparation and standard/QC concentrations.- Include an internal standard with known concentration during sample preparation.- Incorrect sample diluent used.- Check injection volume.
Loss of column performance	<ul style="list-style-type: none">- Check the peak widths and resolution. Test the performance of the column using your standard test for loss in performance.- Replace analytical column.
LC leaks	<ul style="list-style-type: none">- Check for loose fittings post injector on the system.
Mobile phase	<ul style="list-style-type: none">- Check concentration of additives.- If suppressing mobile phase components were used in previous method, clean MS source and flush out LC system.




Undesired Fragmentation

Causes	Solutions
Ion source setting too harsh	<ul style="list-style-type: none">- Check source temperatures are appropriate for analyte.- Check ionisation voltage is appropriate for analyte.
Collision energy too high / low	<ul style="list-style-type: none">- Check and optimise collision cell gas pressure and collision energy.
Ion optics	<ul style="list-style-type: none">- Ensure correct voltage is applied to desolvation line and QArray.



No Peaks

Causes	Solution
MS setting issue	<ul style="list-style-type: none">- Check the method uses appropriate MS settings for the compounds of interest.- Ensure probe is at correct distance from orifice.- Check protrusion of capillary from probe.- Check spray from capillary.- Check ion source temperatures / gas flows are as expected and stable.- Check collision gas pressure is correct and stable.- Ensure analyte elutes within event window set in method.- Check for adducts such as Na⁺, NH₄⁺, K⁺. Consider different glassware, make up new mobile phases.- Check for matrix effects such as ion suppression. Consider different LC methods and sample preparation procedures.- Concentration injected is below limit of detection.- Perform and check MS tune.
LC	<ul style="list-style-type: none">- Check LC outlet tubing is connected to ion source.- No mobile phase flow, possibly purge valve left open.- Purge the system including injector to remove possible air bubbles in pump.- Purge LC system using isopropanol to ensure check valves are working correctly.- Check for crimped or damaged tubing.
Compounds not retained or retained longer than run time in method conditions	<ul style="list-style-type: none">- Check mobile phase composition is correct.- Check correct analytical column type is being used.- Increase run time.- Increase solvent strength.- Check correct flow rate is being achieved.
Sample issues	<ul style="list-style-type: none">- Prepare fresh samples.- Review injection volume in sequence / method.- Ensure the sample is in the correct position in the autosampler.- Check for sample adsorption issue.- Check for air pockets trapped in bottom of vial or well.
Sample flowing to waste	<ul style="list-style-type: none">- Check divert valve settings if applicable.



Basic Steps

Follow these steps to isolate where the problems is.

Check the obvious explanations first and change only one thing at a time!



Check the Basics:

- Power & electrical connections
- Communication cables
- Sample preparation
- Analytical conditions
- Mobile phase preparation
- Needle rinse & seal washes
- Solvent flowing / no air bubbles
- LC pump pressures
- Ion source maintenance
- Roughing pump (oil level & gas ballast)
- MS vacuum
- Argon gas cylinder (level and pressure)
- Gas generator (pressure readings)

Identify the Cause:

- Clearly define the problem
- Review sample and maintenance logs to identify trends in the data or possible problem indicators
- Use a logical sequence of events to isolate possible causes

Document Everything:


- Document typical values such as analyte retention times and normal initial LC operating pressures. Use as a benchmark to indicate deterioration in system performance
- Document all troubleshooting steps and results. This may help you identify and solve the next problem faster
- Always inject a known sample and compare to previous data as a reference to ensure restored performance

Still having problems with your instrument? Let us know at

lc-support@shimadzu.nl


Low Pressure

Causes	Solutions
Partial leak in system	<ul style="list-style-type: none">- Check all connections and retighten any which have leaks.
Flow rate	<ul style="list-style-type: none">- Check the method has the correct flow rate.- Test the flow rate accuracy using a calibrated flow rate meter or collect a specific volume and monitor the time required.- Replace worn out or damaged pump seals.
Method	<ul style="list-style-type: none">- Check if method is using correct temperature and correct solvents.- If a column section valve is used, check correct column selected.
Incorrect column	<ul style="list-style-type: none">- Use correct column with correct dimensions and particle geometry.
Column temperature too high	<ul style="list-style-type: none">- Set adequate column temperature and check no column damage if exceeded column temperature limit.
Airlock in LC tubing	<ul style="list-style-type: none">- Remove tubing from degasser and ensure flow under gravity. Reconnect and purge pumps in isopropanol.
Stuck check valve	<ul style="list-style-type: none">- Purge LC system using isopropanol and ensure check valves are working correctly.- Sonicate the check valves in isopropanol.
Sensor malfunction	<ul style="list-style-type: none">- Repair or replace pressure sensor.




Poor Mass Accuracy (HRAM Instruments)

Causes	Solutions
MS out of tune	<ul style="list-style-type: none">- Perform and check system tune.
TOF Calibration	<ul style="list-style-type: none">- Perform TOF calibration.
Calibration performed incorrectly	<ul style="list-style-type: none">- Ensure sample analytes are within calibration range and adjust if required.
Detector saturation	<ul style="list-style-type: none">- Dilute sample or adjust injection volume.



Changes in MS Resolution

Causes	Solutions
MS out of tune	<ul style="list-style-type: none">- Perform and check tune.



Carry Over

Causes	Solutions
Inappropriate wash settings	<ul style="list-style-type: none">- Check wash solution and wash settings.
Sample concentration	<ul style="list-style-type: none">- Use lower concentrated sample or inject less.
Column contamination	<ul style="list-style-type: none">- Flush column; replace guard columns and analytical column if applicable.
Injector issue	<ul style="list-style-type: none">- Changes in dispensing volume of injector, use a system suitability sample to determine volume changes.- Check batch / method details to ensure the correct volume was programmed.- Increase needle and loop flushing protocols to ensure no carryover from injection.- Purging the injector metering pump.
LC Gradient	<ul style="list-style-type: none">- Insufficient time at strong solvent conditions during gradient program. Increase based on column dimensions.



High Pressure

Causes	Solutions
Flow rate set too high	<ul style="list-style-type: none">- Reduce flow rate setting.
Blocked column	<ul style="list-style-type: none">- Backflush column (if permitted) or replace column.
Incompatible mobile phase (precipitated buffer or immiscible)	<ul style="list-style-type: none">- Use correct mobile phase; wash column and re-equilibrate.
Improper column	<ul style="list-style-type: none">- Use correct column with correct dimensions and particle geometry.
Injector blockage	<ul style="list-style-type: none">- Clear blockage (review needle, loop, valve assembly and HPV outlet).
Guard column / cartridge blockage	<ul style="list-style-type: none">- Replace or remove guard column.
Column in-line filter blockage	<ul style="list-style-type: none">- Replace or remove in-line filter.
Column temperature too low	<ul style="list-style-type: none">- Set adequate column temperature.
Sensor malfunction	<ul style="list-style-type: none">- Repair or replace pressure sensor.
Pump in-line filter blockage	<ul style="list-style-type: none">- Replace in-line filter.
Blocked tubing	<ul style="list-style-type: none">- Replace blocked tubing as necessary.



SHIMADZU
Excellence in Science