


Specification Sheet

Rack Changer Rack Changer /c



	Rack Changer	Rack Changer /c
Sample cooling	No	Yes (see below for details)
Number of samples processed	1152 (12 microtiter plates or deep well plates with 96 wells each)	
Operating conditions	4 - 35°C, 20-85% humidity	
Power requirements	50 VA, 50/60 Hz	350 VA, 50/60 Hz
Dimensions	W 425 x H 420 x D 500 mm	
Weight:	25 kg	32 kg

Rack Changer/c sample cooler specifications

System	Direct cooling system (environment conditions: room temperature below 30°C with humidity less than 70% when the cooler temperature is set to 4°C), dehumidification system built-in
Temperature setting range	4 – 40°C (can be cooled down to 4°C when the room temperature is below 30°C and the humidity is below 70%)
Temperature accuracy	± 6°C